

Oosafe®

IVF LABORATORY DISINFECTANTS

Manufactured in accordance  
with EN ISO 13485:2016


DISIN  
FECT  
ANTS



\*All enveloped viruses from Annex A, EN14476:2013+A1:2015 including HBV, HCV, HDV, HIV, Coronavirus, Influenza, Rabies, Measles virus.

# Oosafe® IVF LABORATORY DISINFECTANTS

## Oosafe® DISINFECTANT FOR CO<sub>2</sub> INCUBATORS AND LAMINAR FLOW HOODS



### High quality

- Ready-to-use solution
- Anticorrosive
- Securely sealed
- Valid for 3 years, even after opening
- Long term disinfection/low evaporation rate
- Fast disinfection effect
- Transparent and odorless

### Health and safety

- MEA and HSSA tested
- Alcohol-free
- No VOC release

### Excellent material compatibility

- Stainless steel
- Plastic
- Glass
- Acrylic

### Application:

Follow laboratory's surface cleaning procedure before applying the disinfectant.

**Spray:** Spray and distribute evenly on the surface with a sterile cloth. Wait 15 minutes for the disinfection effect. Wipe the surface to remove any residue.

**Wipes:** Wipe the surface evenly and wait 15 minutes for the disinfection effect. Close the lid after use. Do not reuse wipes.

### Intended use:

Disinfection of CO<sub>2</sub> incubators and laminar flow hoods.



Product reference:  
**OODIH-0050**

Oosafe® Disinfectant for CO<sub>2</sub> Incubators and Laminar Flow Hoods, 50 ml in spray bottle



Product reference:  
**OODIH-1000**

Oosafe® Disinfectant for CO<sub>2</sub> Incubators and Laminar Flow Hoods, 1 Liter bottle with spray



Product reference:  
**OODIH-5000**

Oosafe® Disinfectant for CO<sub>2</sub> Incubators and Laminar Flow Hoods, 5 Liters refill container



Product reference:  
**OODW-70**

Oosafe® Disinfection Wipes, 70 wipes per bucket

# Oosafe® IVF LABORATORY DISINFECTANTS

## Oosafe® LABORATORY SURFACE DISINFECTANT

### High quality

- Securely sealed
- Valid for 3 years, even after opening
- Long term disinfection

### Health and safety

- MEA and HSSA tested
- Clear - slightly yellow
- Alcohol-free
- No VOC release

### Excellent material compatibility

- Glass, Plastic, Rubber, Stainless steel, Ceramic, Linoleum, Acrylic, Cotton

### Features

- Contains non-ionic surfactant to clean and quaternary ammonium compounds to disinfect surfaces at the same time.

### Application:

Dilute disinfectant with water according to the table recommendations and apply the diluted cleaning agent on laboratory floors and surfaces. Wait after application for the required disinfection effect and let the surface dry out. If used for laboratory floors, please rinse otherwise it might be slippery.

### Intended use:

For disinfection and washing of laboratory surfaces and floors.



Product reference:  
**OODSF-00050**

Oosafe® Laboratory Surface Disinfectant,  
50 ml in bottle



Product reference:  
**OODSF-02000**


Oosafe® Laboratory Surface Disinfectant,  
2 Liters bottle with measuring cup



Product reference:  
**OODSF-10000**

Oosafe® Laboratory Surface Disinfectant,  
10 Liters refill container

\*All enveloped viruses from Annex A, EN14476:2013+A1:2015 including HBV, HCV, HDV, HIV, Coronavirus, Influenza, Rabies, Measles virus.






\*All enveloped viruses from Annex A, EN14476:2013+A1:2015 including HBV, HCV, HDV, HIV, Coronavirus, Influenza, Rabies, Measles virus.

# Oosafe® IVF LABORATORY DISINFECTANTS

## Oosafe® HAND DISINFECTANTS



Product reference:  
**OODH-0050**

Oosafe® Hand Disinfectant,  
50 ml in spray bottle



Product reference:  
**OODH-0500**

Oosafe® Hand Disinfectant,  
500 ml bottle with pump



Product reference:  
**OODH-5000**

Oosafe® Hand Disinfectant,  
5 Liters refill container

### High quality

- Ready-to-use solution
- Securely sealed
- Valid for 3 years, even after opening
- Long term disinfection/low evaporation rate
- Fast disinfection effect
- Transparent and odorless

### Health and safety

- Non-sensitizing
- Non-irritating
- MEA and HSSA tested
- Alcohol-free
- No VOC release
- Keeps your skin soft

### Application:

Follow laboratory's hygienic procedure for washing hands before using hand disinfectant. Apply minimum 3 ml of disinfectant and rub hands until they are completely dry.

### Intended use:

For disinfection of hands.

#### ORDER INFORMATION

##### Oosafe® Disinfectant for CO<sub>2</sub> Incubators and Laminar Flow Hoods

Product reference	Product name
OODIH-0050	50 ml in spray bottle
OODIH-1000	1 Liter bottle with spray
OODIH-5000	5 Liters refill container
OODW-70	70 wipes per bucket

##### Oosafe® Laboratory Surface Disinfectant

Product reference	Product name
OODSF-00050	50 ml in spray bottle
OODSF-02000	2 Liters bottle with measuring cup
OODSF-10000	10 Liters refill container

##### Oosafe® Hand Disinfectant

Product reference	Product name
OODH-0050	50 ml in spray bottle
OODH-0500	500 ml bottle with pump
OODH-5000	5 Liters refill container

f facebook.com/SparMED

t twitter.com/SparMED

in linkedin.com/company/sparMED-ApS

**sparmed**

SparMED ApS  
Ryttermarken 2, 3520 Farum, DENMARK  
Phone number: +45 39-40-25-03  
MK\_4.3.B10 rev. 20.06.2023

