



Instruction for Use

Oosafe® Disinfectant

DocNo. OODIS-21 Rev. No. 07 Rev. Date 31.05.2022



Scope of Application

The product Oosafe® Disinfectant is a surface disinfectant based on quaternary ammonium compounds suitable for medical devices - dental and surgical instruments. The product is intended to disinfect non-porous surfaces and is effective against bacteria, yeasts, viruses, and mycobacteria (see table below). Due to no Volatile Organic Compound (VOC) release and no fragrance it is a favourable disinfectant for medical devices in laboratories, including those with a low ventilation level.

Disinfection of a surface is carried out to prevent the spread of infectious diseases and it is used to inhibit the activity or destroy harmful microorganisms and to protect people during the pre-treatment of medical instruments. The surface should be clean before disinfection can occur.

Application

- **Recommended waiting time:** 15 min (see table)
- **Recommended use:**
- **Apply undiluted**
- **A method - Spray:** Make sure the surface to be disinfected is clean. Spray and distribute the disinfectant evenly on the surface with a sterile cloth. Wait 15 min for the disinfection effect.
- **B method - Wipes:** Make sure the surface to be disinfected is clean. Wipe the surface evenly and wait 15 min for the disinfection effect.

Properties

- Transparent and odourless
- Ready-to-use
- Alcohol-free
- Non-VOC releasing
- Effective against mycobacteria
- Effective against enveloped viruses including HIV, Hepatitis B and C, Herpes simplex, Influenza A, Ebola, also Adenovirus and Poliovirus
- Does not cause corrosion of metals

Durability

After opening Oosafe® Disinfectant bottles or containers, mixture maintains stability for 3 years under the proper storage conditions.

Storage

- Store products in a dry place (ideal humidity range of 30 – 50%) at a temperature between 15°C to 30°C
- Avoid direct sunlight

Efficacy

*Efficacy is tested for clean conditions

**Effectiveness against all enveloped viruses including HIV, Hepatitis B and C, Herpes simplex, Influenza A, Ebola

Activity Spectrum	Method*	Contact Time
Bactericidal <i>Pseudomonas aeruginosa, Staphylococcus aureus, Enterococcus hirae</i>	EN 13727+A1 EN 13727+A2 EN 1040 EN 14561	5 min 5 min 15 min 15 min
Yeasticidal <i>Candida albicans</i>	EN 1275 EN 13624 EN 14562	5 min 15 min 15min
Mycobactericidal and Tuberculocidal <i>Mycobacterium avium, Mycobacterium terrae</i>	EN 14348 EN 14563	15 min 15 min
Virucidal <i>Vaccinia virus**</i>	EN 14476+A1	15 min
<i>Poliovirus</i>	EN 14476	15 min
<i>Adenovirus</i>	EN14476+A1	5 min

Packaging Type

Packaging Type	Product code
1. 50 ml, spray bottle	OOMD-00050
2. 250 ml, spray bottle	OOMD-00250
3. 1L, spray bottle: individually packed in the carton	OOMD-01000
4. 5L, refill container, individually packed in the cardboard protection	OOMD-05000
5. 10, refill container, individually packed in cardboard protection	OOMD-10000
6. Wipes, bucket with 70 wipes	OOMDW-120

Safety Information

Dangerous reaction:

No dangerous reactions known if used and stored according to specifications

Toxicological information:

- On the skin:
No irritating effects known.
- In the eye:
Rinse immediately with water if irritation arises. Consult the doctor if irritation persists.

Doc. No. OODIS-21 Rev. No. 07 Rev. Date 31.05.2022









- Ingestion:
Rinse mouth and drink plenty of water. Consult a doctor.
- Inhalation:
No special measures required. Supply fresh air and consult a doctor if any discomfort arises.
- Sensitization:
No sensitizing effects known.
- Disposal considerations:
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

For additional information about industrial safety refer to Oosafe® Disinfectant safety data sheet.

Notice to User:

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user is established.

Symbols on the product's primary packing label

Symbol	Description of the symbol
1. 	Date of manufacture
2. 	Keep dry
3. 	Consult instructions for use or consult electronic instructions for use
4. 	Manufacturer
5. 	Use by date
6. 	Temperature limit
7. 	Batch code
8. 	Conformity Mark with identification number of the notified body