


Oosafe[®]
CRYOWARE

CRYO





Oosafe® CRYOWARE



OOSAFE® COLOR CANES

Oosafe® Color Canes are designed for safe storage of cryo vials and straws (inside goblet/visiotube). Unlimited and unique ID combinations can be made for the stored samples together with different colors of goblets and printed/non-printed tags. This great feature helps to differentiate each samples in LN₂ tanks and minimizes the risk of mix-ups.

The whole range of canes is made of aluminum. There are different types of canes to be used with different sizes of goblets and 1-2 ml cryo vials.



ORDER REFERENCE

All color canes are available for sale only by pre-order of the minimum quantity
1.000 pcs = 20 packs of 50.

	PRODUCT	SILVER	BLUE	GOLD	GREEN	PINK	PURPLE	RED
1	Oosafe® Color Cane for 1-2 ml vials with 5 security shelves, pack of 50 pcs. Holds 5 pcs of 1-2 ml cryo vials.	OOCW-CC01	OOCW-CC03	OOCW-CC04	OOCW-CC05	OOCW-CC06	OOCW-CC07	OOCW-CC08
2	Oosafe® Color Cane for 13 mm diameter goblet (visiotube) with 1 security shelf, pack of 50 pcs. Holds 1 pc of 13 mm diameter goblet (visiotube).	OOCW-CC11	OOCW-CC13	OOCW-CC14	OOCW-CC15	OOCW-CC16	OOCW-CC17	OOCW-CC18
3	Oosafe® Color Cane for 12 mm diameter goblet (visiotube) with 1 security shelf, pack of 50 pcs. Holds 1 pc of 12 mm diameter goblet (visiotube) or 6 pcs of 1-2 ml cryo vials.	OOCW-CC21	OOCW-CC23	OOCW-CC24	OOCW-CC25	OOCW-CC26	OOCW-CC27	OOCW-CC28




OOSAFE® PRINTED & NON-PRINTED TAGS FOR CANES


The whole range of colored printed & non-printed ID tags are made of aluminum and the color line is powder coated. They are safe to use in the storage containers in the temperature down to -100°C. Prints on tags are durable and permanent, unlike handwritten tags with permanent cryo marker pens. This feature makes printed tags the safest choice for storage inside LN₂ tanks.



Oosafe® Non-printed tags,
pack of 50 pcs



Oosafe® Printed tags,
pack of 35 pcs *



Standard printed tags are numbered from 1 to 35 but printing can also be customized.
*Please contact SparMED regarding the minimum purchase order quantity required for customized printed tags.

ORDER REFERENCE

All tags for canes are available for sale only by pre-order of the minimum quantity
1.000 pcs = 20 packs of 50.

PRODUCT	YELLOW	ORANGE	GREEN	WHITE	PINK	BLUE	PURPLE	SILVER	RED	BROWN
Oosafe® Non-printed tags, pack of 50 pcs	OOCW-PT01	OOCW-PT02	OOCW-PT03	OOCW-PT04	OOCW-PT05	OOCW-PT06	OOCW-PT07	OOCW-PT08	OOCW-PT09	OOCW-PT10
Oosafe® Printed tags, pack of 35 pcs	OOCW-PT11	OOCW-PT12	OOCW-PT13	OOCW-PT14	OOCW-PT15	OOCW-PT16	OOCW-PT17	OOCW-PT18	OOCW-PT19	OOCW-PT20







OOSAFE® SLEEVES FOR CANES

The unique transparent color design provides security and visibility while handling and storing cryo vials and goblets (visotubes) inside LN₂ tanks.

Different colors allow the user to separate or group patients' samples, as well as to locate them faster inside LN₂ tanks.

Slightly smaller diameter, compared to existing sleeves on the market, allows usage of cryo canes with both 12 mm and 13 mm diameter goblets (visotubes) and take up less space inside LN₂ tank canisters.



ORDER REFERENCE

PRODUCT	Transparent	Transparent Blue	Transparent Green	Transparent Yellow	Transparent Pink	Transparent Purple
Oosafe® Sleeve for Canes	OOCW-CS01	OOCW-CS02	OOCW-CS03	OOCW-CS04	OOCW-CS05	OOCW-CS06



facebook.com/SparMED



twitter.com/SparMED



linkedin.com/company/SparMED-ApS

sparmed

SparMED ApS

Ryttermarken 2, 3520 Farum, DENMARK

Phone number: +45 39-40-25-03

MK_4.3.B4 rev. 20.06.2023

