

Oosafe® IVD Medical Device **PLASTIC WARE**



Manufactured in ISO 9001:2015
and EN ISO 13485:2016 certified company.

Oosafe® IVD Medical Device PLASTIC WARE



TOP QUALITY IVD MEDICAL DEVICE PLASTICWARE APPROVED BY SCIENTISTS AND SUPPORTED BY THE TEST RESULTS FOR EACH LOT


- Tested for each LOT:
 - One cell stage Mouse Embryo Assay (MEA > 80% Blastocyst)
 - Non-pyrogenic, Endotoxin (LAL < 0.03 EU/device)
 - Human Sperm Survival and Motility (SMI ≥ 0.75)
- Surface treated options for a consistent hydrophilic surface
- Packaged in peel or tear-open packaging for user convenience
- Easy Grip and Label Area Design
- Sterility assurance level (Irradiation) 10^{-6}
- CE marked according to 98/79/EC Directive on *In Vitro* Diagnostic Medical Devices
- Made of tested Raw Materials
 - Polystyrene (PS)
 - Polypropylene (PP)
- Vacuum packed to prevent scratches during transportation



Oosafe® 4 Well Dish TREATED SURFACE

MULTIPACK PS

4 pcs/pack
120 pcs/case



Product reference:
OOPW-FW03

Cut corner design ensures proper orientation
and prevents mix-ups while handling
Well volume 1.8 mL
Volume between wells 7 mL


BREATHABLE PACKAGING - NO OFF GASSING NEEDED



Oosafe® 4 Well Dish NON-TREATED SURFACE

MULTIPACK PS

4 pcs/pack
120 pcs/case



Product reference:
OOPW-FW04

Cut corner design ensures proper orientation
and prevents mix-ups while handling
Well volume 1.8 mL
Volume between wells 7 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED



Oosafe® 6 Well Dish with Straw Holder NON-TREATED SURFACE

SINGLEPACK PS

1 pc/pack
240 pcs/case



Product reference:
OOPW-SW02

Conical flat bottom. Hands-free design: the tip of the straw
and the bottom of the well are at the same focus level
Well volume 0.4 mL
Volume between wells 3 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED



Oosafe® 6 Well Dish with Straw Holder NON-TREATED SURFACE

MULTIPACK PS

4 pcs/pack
120 pcs/case



Product reference:
OOPW-SW03

Conical flat bottom. Hands-free design: the tip of the straw
and the bottom of the well are at the same focus level
Well volume 0.4 mL
Volume between wells 3 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® Center Well Dish with 2 Compartments

ALL AROUND GRIP
NON-TREATED SURFACE

MULTIPACK PS

10 pcs/pack
500 pcs/case



Product reference:
OOPW-CW03

Multi-purpose dish for cell handling applications
Center well volume 3 mL
Each compartment volume 5 mL

VACUUM PACKED TO PREVENT SCRATCHES

Oosafe® Center Well Dish with 2 Compartments

NON-TREATED SURFACE

MULTIPACK PS

10 pcs/pack
500 pcs/case



Product reference:
OOPW-CW05

Multi-purpose dish for cell handling applications
Center well volume 3 mL
Each compartment volume 5 mL


VACUUM PACKED TO PREVENT SCRATCHES

Oosafe® Center Well Dish

NON-TREATED SURFACE

MULTIPACK PS

10 pcs/pack
500 pcs/case



Product reference:
OOPW-CW08

New design with one outer chamber
Center well volume 3 mL
Outer compartment volume 10 mL

VACUUM PACKED TO PREVENT SCRATCHES

Oosafe® 50 mm Dish, Thin Bottom

NON-TREATED SURFACE

MULTIPACK PS

10 pcs/pack
500 pcs/case



Product reference:
OOPW-IC06

Extra thin observation area for better image quality
Thickness of observation area: 0.9 mm
Thickness of outer ring: 1 mm
Volume 12 mL


VACUUM PACKED TO PREVENT SCRATCHES

Oosafe® 50 mm Dish, Thin Bottom

NON-TREATED SURFACE

MULTIPACK PS

10 pcs/pack
500 pcs/case



Product reference:
OOPW-IC07

Design with two channels for sperm selection.
Channel volume: 15.8 µL for each
Volume 13.3 mL

VACUUM PACKED TO PREVENT SCRATCHES

Oosafe® 35 mm Dish, High Wall

TREATED SURFACE

MULTIPACK PS

10 pcs/pack
500 pcs/case



Product reference:
OOPW-TF05

Multi-purpose dish for cell culture and microscopy
applications, high wall with enough space
for labels between the grips
Volume 8 mL

VACUUM PACKED TO PREVENT SCRATCHES

Oosafe® 35 mm Dish

NON-TREATED SURFACE

MULTIPACK PS

10 pcs/pack
500 pcs/case



Product reference:
OOPW-TF06

Multi-purpose dish for cell culture and microscopy
applications, high wall with enough space
for labels between the grips
Volume 8 mL

VACUUM PACKED TO PREVENT SCRATCHES

Oosafe® 60 mm Dish,

NON-TREATED SURFACE

MULTIPACK PS

10 pcs/pack
500 pcs/case



Product reference:
OOPW-ST06

Multi-purpose dish for cell culture
and microscopy applications
Volume 25 mL

VACUUM PACKED TO PREVENT SCRATCHES



Oosafe® 100 mm Dish
NON-TREATED SURFACE

MULTIPACK PS

10 pcs/pack
250 pcs/case

Product reference:
OOPW-HD10

Multi-purpose dish for cell culture
and microscopy applications
Volume 105 mL

VACUUM PACKED TO PREVENT SCRATCHES



Oosafe® Synergy Dish
TREATED SURFACE

MULTIPACK PS

4 pcs/pack
120 pcs/case



Product reference:
OOPW-SD04

Unique design with micro-wells for group culture
Volume of wash well 56 µl
Volume of group culture well 160 µl
Total volume 15 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® Sperm Collection Cup, 80 mL

NON-TREATED SURFACE

SINGLEPACK PP

1 pc/pack
70 pcs/case



Product reference:
OOPW-SC01

Wide opening for easy sampling
Minimum volume measured 2 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® Andrology Tube, 5 mL

NON-TREATED SURFACE

MULTIPACK PS

10 pcs/pack
500 pcs/case



Product reference:
OOPW-AT10

Push-on cap with dual position to enable venting
Minimum volume measured 0.5 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® Snap-Cap Tube, 14 mL

NON-TREATED SURFACE

MULTIPACK PS

10 pcs/pack
500 pcs/case



Product reference:
OOPW-OT10

Push-on cap with dual position to enable venting
Minimum volume measured 0.1 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® PS Centrifuge Tube, 15 mL

NON-TREATED SURFACE

MULTIPACK PS

5 pcs/pack
500 pcs/case



Product reference:
OOPW-CT05

Suited for procedures requiring high optical clarity
Minimum volume measured 0.1 mL

BREATHABLE PACKAGING - NO OFF GASSING NEEDED

Oosafe® IVD Medical Device PLASTIC WARE



| Product reference | Product name | Surface treatment | Packaging |
|-------------------|---------------------------------------------------------------|-------------------|------------|
| OOPW-FW03 | Oosafe® 4 Well Dish | Treated | MULTIPACK |
| OOPW-FW04 | Oosafe® 4 Well Dish | Non-treated | MULTIPACK |
| OOPW-SW02 | Oosafe® 6 Well Dish with Straw Holder | Non-treated | SINGLEPACK |
| OOPW-SW03 | Oosafe® 6 Well Dish with Straw Holder | Non-treated | MULTIPACK |
| OOPW-CW03 | Oosafe® Center Well Dish with 2 Compartments | Non-treated | MULTIPACK |
| OOPW-CW05 | Oosafe® Center Well Dish with 2 Compartments, Label Area Grip | Non-treated | MULTIPACK |
| OOPW-CW08 | Oosafe® Center Well Dish | Non-treated | MULTIPACK |
| OOPW-IC06 | Oosafe® 50 mm Dish, Thin Bottom | Non-treated | MULTIPACK |
| OOPW-IC07 | Oosafe® 50 mm Dish, Thin Bottom | Non-treated | MULTIPACK |
| OOPW-TF05 | Oosafe® 35 mm Dish, High Wall | Treated | MULTIPACK |
| OOPW-TF06 | Oosafe® 35 mm Dish | Non-treated | MULTIPACK |
| OOPW-ST06 | Oosafe® 60 mm Dish | Non-treated | MULTIPACK |
| OOPW-HD10 | Oosafe® 100 mm Dish | Non-treated | MULTIPACK |
| OOPW-SD04 | Oosafe® Synergy Dish | Treated | MULTIPACK |
| OOPW-AT10 | Oosafe® Andrology Tube, 5mL | Non-treated | MULTIPACK |
| OOPW-OT10 | Oosafe® Snap-Cap Tube, 14 mL | Non-treated | MULTIPACK |
| OOPW-CT05 | Oosafe® PS Centrifuge Tube, 15 mL | Non-treated | MULTIPACK |
| OOPW-SC01 | Oosafe® Sperm Collection Cup, 80 mL | Non-treated | SINGLEPACK |

All plasticware products are packed in a secure white cardboard box.

Oosafe® IVD Medical Device PLASTIC WARE

sparmed



SparMED ApS
Ryttermarken 2, 3520 Farum, DENMARK
Phone number: +45 39-40-25-03



MK_4.3.B2 rev. 06.10.2023

